

ORTHOLESE



Animal Orthotics

Wearing the orthosis

Your Veterinarian or healthcare professional will show you how to put your dog's orthosis on. Make sure that it is fitting correctly and that you fasten the straps firmly on so that the orthosis is secure.

Your Veterinarian or healthcare professional will tell you how long the orthosis should be worn for and provide you with a wearing regime. Make sure you follow the regime at all times as continuity of wear is critical to the effectiveness of the treatment.

It can take time for your dog to get used to wearing their orthosis. Whilst your veterinarian or healthcare professional will ensure that it is as safe and comfortable as possible to wear, it may slightly alter the way they move. This is completely normal, however we do advise you gradually build up the length of time their orthosis is worn to help them become accustomed to it, and to ensure it will not cause ill effects such as skin irritation or damage.

IMPORTANT NOTE REGARDING THE INFORMATION PROVIDED

If you see or suspect any new or prolonged symptoms in your pet, please contact your Vet/Physio/healthcare professional for a full examination/diagnosis. Please do not carry out any treatments without prescriptions from the experts.